

## FOR RELIABLE, RUGGED AND PRECISE BALANCES, YOU CAN PUNCH ON THE PAN, NOT POSSIBLE TO DESTROY THE INTERNAL WEIGHING SYSTEM




## FOR RELIABLE, RUGGED AND PRECISE BALANCES

All Gibertini balances are made with one special ANTISHOCK SYSTEM that allows it the use in hard environments, not need specialized staff

Figure 1 shows how the load cell of our balances are equipped with 4 springs, the shock it is absorbed! These springs not present in the common scales, figure 2, in which the plate is mounted directly to the load cell and the shock can cause damage.

### Balance Gibertini



### Other Balances not "Gibertini"

